

# Dr's Cam Air

*Enjoy Wireless*

Wireless AutoFocus Intraoral Camera with Mouse



















Good Doctors Wireless Solution	Caries Detector	Wireless Intraoral Camera	Wired Intraoral Camera		Apex Locator	Curing Light	
IC-AD	IC-WHCD100 IC-DCDU75	IC-WHT60 <b>IC-DCAM70</b> IC-EZAM80	IC-WHMA90 IC-WHW61 IC-DCAU71	IC-EZAU81 IC-DCCU74 IC-EZBU82	AL-DFA20 AL-DF10	CL-DP40 CL-DP41 CL-DL10	CL-DC20

# Enjoy Wireless

## Wireless AutoFocus Intraoral Camera with Mouse

### 1 Features of Dr's cam Air

-  **High speed digital wireless**  
High speed digital wireless provide more clear live image without interruption.
-  **Macro to Full arch**  
Camera is able to have Macro to Full arch Image.
-  **Auto Focus**  
Real Autofocus function shows clear images.
-  **FHD image support**  
The device shows FHD (1080p) in live image supported by a 5M sensor.
-  **A Simple selection and Easy connection**  
Select the Connection signal on the camera.
-  **Infinite Channel**  
Unlimited Channel Support.
-  **Real Plug and Play**  
Not required Driver Software.
-  **Mac Support**  
Mac and PC support with wireless.
-  **Wireless mouse**  
Built in Gyro mouse.
-  **OLED screen**  
Bright and Clear OLED screen indicates various device status.
-  **280° twist head**  
Easy to get approach anywhere in oral cavity.
-  **Adjustable LED brightness**  
Can adjust LED brightness to avoid light reflecting.
-  **3 Color Impressions**  
Possible to Select Color Impression(Original, Warm, Cool)
-  **Software compatibility**  
Able to be used with almost any dental imaging software in the market all over the world.
-  **Twain driver**  
Twain Driver Provided.
-  **Big capacity battery**  
Long time use with big capacity battery.



### 2 Specification

Model name	IC-DCAM70	
Camera	Operation voltage	3.7V
	Size	Approx. 255 x 30 x 28 mm
	Weight	Approx. 89g
	Resolution	FHD (1920x1080)
	Sensor	5M Sensor CMOS / Liquid
	Charger	Operation voltage
Charger	Size	Approx. 110 x 105 x 68 mm
	Weight	Approx. 100g
	Disposable cover	Size
Disposable cover	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

### 3 Minimum System Requirement

PC system:	3.3 Ghz CPU, 4GB RAM
OS:	Windows 7
USB port:	USB 2.0

