

# Dr's Light Clever

## Incomparable Dual Wave Light Quality

Dr's Light Clever is a new conceptual curing light applied many advantages for user's demands.  
Just press the button and experience new light technology with Dr's Light Clever.













# Incomparable Dual Wave Light Quality


Dr's Light Clever is a new conceptual curing light applied many advantages for user's demands.  
Just press the button and experience new light technology with Dr's Light Clever.



## 1 Features of Dr's Light Clever Dual

-  **Dual wavelength LED**  
Usable all kinds of materials in the field with 405nm & 460nm peak wavelength LED.
-  **Focused light**  
Increase efficiency of light.
-  **Autoclaveable and interchangeable Guide-cap**  
Guide-cap is autoclaveable and interchangeable to prevent cross infection.
-  **4 curing light modes**  
User can select High, Low, Softstart, Pulse mode as user wants.
-  **Time setting available (Softstart is fixed for 20 sec.)**  
User can select time setting as needs.
-  **330° head twist and slim head**  
It makes easy to get approach into the oral.
-  **Angled handpiece design**  
Most familiar design for user.
-  **OLED screen**  
It shows device states on the OLED screen with clearly.
-  **Interchangeable battery pack**  
User can interchange the battery easy.  
(Doesn't need to send device to factory to change the battery)
-  **Big capacity battery (Samsung cell)**  
Big capacity battery makes long time use.

## 2 Comparison with Other products




## 4 Specification


Model : CL-DC20 / Wave length : 400~490(nm) / Guide-cap : Autoclaveable /  
Power source : Lithium-ion-battery / Light source : High Power Dual wavelength LED /  
Weight (Handpiece) : Approx. 130g

## 3 Various curing modes


Selectable light intensity and time as you want



## 5 Accessory



Guide-cap (x 3 pcs)



Battery